



Quick Guide



01 Mounting a sample chamber

- A The sample chamber is removed from the storage tubes. Note that this should first be done inside the test section.
- B Insert the sample chamber in the top of the AeroCollect® as shown.
- C By pressing the green push-button the sampling is initiated. The red LED will light up as the 400 V voltage drop is ramping up. Once the green LED turns on the sampling procedure has started. Note that the instrument will not start sampling if a sample chamber is not properly inserted (In that case the red LED will remain lighted and the green will not turn on).



02 Sample collection procedure

When the sample chamber has been correctly inserted the collection procedure can be initiated.

The bacteria and virus are collected from the air inside the pig facilities. The sample should be collected from the entire section by walking along the central walkway with the AeroCollect® Instrument held over the pens.

Sampling should be conducted for a period of minimum 15 minutes while at the same time ensuring that the sample are collected from all parts of the section of interest.



03 Stopping the collection procedure

- A After covering the entire house the sampling procedure is terminated by pressing the red push-button on the front of the sampler.
- B When the AeroCollect® has been stopped the sample chambers can be removed by lifting the tip of the sample chamber while pulling it out.
- C The filled collection chamber is reinserted into the storage tube. Please ensure to carefully label the storage tube with sample section identification for future reference.

04 The collected samples are shipped for analysis to your test service laboratory

Questions and further comments are welcome to
FORCE Technology
AeroCollect
Venlighedsvej 4
2970 Hørsholm
Denmark

Tel. +45 43 25 14 00
aerocollect@force.dk
aerocollect.dk